

LABScreen™ 2.3

Advance Preview


Graphical Display Module

Several new graphic features have been added to the software supported in the 2.3 release. Please take a few moments to review this document to see what has changed. Then, please run the software and let us know what you think of it. Send your comments via e-mail to Nori Sasaki at Nsasaki@OneLambda.com . Please cc: Hillary Janikula at Hjanikula@OneLambda.com.

New Features in This Prototype

- New Options for cutoff adjustment
 - Grid only
 - Grid and graph
 - Graph only

USING THE GRAPHICAL DISPLAY



The Adjust Cutoff functions button has been enhanced so that it now contains three options:

- the "table" on the button on the left indicates a "grid-only" data display.
- The "adjust cutoffs" button represents a "grid and graph" selection.
- The "bar graph" on the right indicates a "graph only" selection.

Using the Grid-only Display

If you choose the “grid-only” button on the Analysis Results screen, this is what appears:

Revise Cutoff/RawData for <90> CatID <53400 001 >

From the Previous Settings										Revised	
Beads	Run	RawData	MBG Ratio	Der. MBG	Prev. MBG	Normalized BG	BG Value	Run	Count	Specificities	
929	2	7423.500	2.334	3165.17	3165.17	3185.17	180.22	2	186	942 B945	
930	1	780.000	0.350	2227.95	2227.95	2227.95	112.07	1	174	944 B944	
931	1	1214.000	0.526	2308.06	2308.06	2380.06	116.14	1	165	945 B945	
932	1	681.000	0.398	2207.08	2207.08	2287.00	111.02	1	190	949 B944	
933	1	1529.000	0.828	2117.42	2117.42	2117.42	106.51	1	111	951 B944	
934	1	1443.000	0.508	2700.01	2700.01	2780.81	139.08	1	75	952 B944	
935	1	1584.000	0.461	2915.68	2915.68	2915.60	146.66	1	190	953 B944	
936	2	8246.000	2.400	2598.51	2598.51	2598.51	130.78	2	121	955 B945	
937	1	2689.000	1.196	2248.03	2248.03	2248.03	112.08	1	140	957 B944	
938	2	9077.000	2.304	2203.78	2203.78	2283.70	110.85	2	111	960 B945	
939	1	989.000	0.374	2322.78	2322.78	2322.78	116.88	1	271	962 B945	
940	1	1471.000	0.538	2621.68	2621.68	2621.68	131.87	1	229	964 B945	
941	1	980.500	0.342	2487.17	2487.17	2487.17	123.60	1	208	972 B945	
942	1	978.000	0.347	2648.43	2648.43	2648.43	133.12	1	234	975 B945	
983	1	670.000	0.383	7342.28	7342.28	7342.28	389.33	1	613	A28	
984	1	2244.000	0.497	4808.38	4808.38	4888.38	241.77	1	729	A08	

Adjust ALL BG Values from Previous Settings [8:30] [Apply] [Defaults] [Update] [Cutoff] [Save] [Print]

NC = 1080.0 PC = 5847.0 PCNC = 2.94

Update BG Values: [Update from] [Save to] [Copy]


Sort Data By: [MBG Ratio] [Reaction] [Bead ID] [RawData]

The functions available on this screen are the same as in the previous release.

Using the Grid and Graph Selection

On the Results screen, click **Adjust Cutoff**.

The Revise Cutoff/RawData screen appears:




This graphical representation offers the following:

- The X axis shows the bead number and specificity for each bead.
- The Y axis shows the ratio to the negative background.
- Y axis scaling will adjust, depending on the high and low values.
- The X axis display can be ordered by a) by bead id, b) by ratio, c) by raw data, d) by reaction. (Specificity may be added later, depending on screen space.)
- You can adjust cutoff values in 1 of 2 ways. You can use the mouse to click and drag the cutoff line; or, you can click **Cutoff** to make numerical entries.
- The graph indicates the definition of positive.

Using the Graph Only Selection

If you select **Graph Only**, the entire layout is devoted to the graphical display of the data.



To adjust cutoffs with your mouse:


1. point to the Cutoff line you want to adjust, and click with the left mouse button.
2. Drag the line to where you want to set the new cutoff, and release the mouse button.

You can also adjust the cutoffs numerically from this screen. Click **Cutoff** and make the adjustments on the Revise Cutoff screen.

Printing the Graph

Printing the graph requires taskbar to be activated.

1. Right-click the graph area.
2. On the floating menu that appears, select **Taskbar**.



3. A printer icon will appear on the task bar. Click it to print the graph.

