


Downloading Templates from Download.onelambda.com

To download individual Templates, click on the links below.

<http://download.onelambda.com/pub/softwareupdate/labscreen/IS22template>
http://download.onelambda.com/pub/software_update/labtype/template_IS22/

Right-click the .idt file name below and choose **Save Target as** or **Save as**.



The browser may default to save it as a WinZip file. Be sure to save it as a .idt file. To do this: open the **Save As Type** list box, and then select **All files**. When you download the file it will be properly formatted as an .idt file. (see Figure 1).




Figure 1: Correcting the **Save As** File type.